

MS6-Pro

LCD Digital Magnetic Stirrer

- High resolution LCD displays actual stirring speed
- Brush less DC motor is maintenance free
- Aluminum covered ceramic coating, corrosion resistant and easy to clean
- Digital speed control up to 1500rpm
- Maximum stirring quantity of H₂O at 5L



Specifications	MS6-Pro
Work plate Dimension [W x D]	140x140mm
Work plate material	Aluminum cover with ceramic coating
Motor type	Brushless DC motor
Stirring positions	1
Max. stirring quantity[H ₂ O]	5L
Max.magnetic bar[length]	40mm
Speed range	200-1500rpm
Speed display	LCD
Protection class	IP21
Motor input Power	2.4W
Motor output Power	1.8W
Power	10W
Voltage,Frequency	100-240V, 50/60Hz
Dimension[WxDxH]	320x180x108mm
Weight	2.2kg
Permissible ambient temperature and humidity	5-40°C, 80%RH