

# MS7-S

## 7" Square Magnetic Stirrer

- Speed range of 0-1500rpm
- Maximum stirring quantity of H<sub>2</sub>O up to 10L
- Ceramic work plate provides excellent corrosion resistance



Specifications	MS7-S
Work plate Dimension [W x D]	184x184mm (7 inch)
Work plate material	Glass ceramic
Motor type	Shaded pole motor
Stirring positions	1
Max. stirring quantity[H <sub>2</sub> O]	10L
Max.magnetic bar[length]	80mm
Speed range	0-1500rpm
Speed display	Scale
Protection class	IP21
Motor input Power	15W
Motor output Power	1.5W
Power	30W
Voltage,Frequency	100-120/200-240V, 50/60Hz
Dimension[WxDxH]	215x360x112mm
Weight	3.8kg
Permissible ambient temperature and humidity	5-40°C, 80%RH

