

# FlatSpin

## Ultra-flat Compact Magnetic Stirrer

With excellent stirring functionality, the compact shape of FlatSpin takes less bench space and makes storage easier.



Specifications	FlatSpin
Work plate Dimension	Ø90mm
Work plate material	PET
Stirring positions	1
Max. stirring quantity[H <sub>2</sub> O]	0.8L
Max.magnetic bar[length]	40mm
Speed range	15-1500rpm IP21
Protection class	3W
Output Power	5W
Power	100-240V,
Voltage,Frequency	50/60Hz
Dimension[WxDxH]	175x120x1mm
Weight	0.3kg
Permissible ambient temperature and humidity	5-40°C,80%RH

## Features

- Speed range 15 to 1500rpm
- Automatic reverse rotation for excellent mixing
- PET top work plate provides excellent chemical resistance
- Motorless technology provides ultra silent and vibration-free operation
- Perfect magnetic stirrer for use during titration (bottle-top burette) in analytical chemistry