

SUPERIOR QUALITY FOR OPTIMUM PERFORMANCE

Accumax Centrifuge Tubes

ACCUMAX

Flawless design, highest standards of execution defines Accumax Centrifuge Tubes. Manufactured in class 100K clean room facility with Swiss-made diamond finish mould, It offers optimum convenience & durability.



Tamper Proof
Graduation Mark



Leak Proof
Screw Cap



Ample Writing Space
On Cap & Tube Surface



Certified DNase, RNase,
Endotoxin, Pyrogen &
Human DNA Free



15 ML
17000 g



50 ML
20000 g



50 ML
Self Standing

For more details visit us at www.caresourceglobal.com

AN ISO 9001 : 2015 & ISO 13485 : 2016 CERTIFIED MFG. COMPANY

*We also execute OEM projects. Specifications can be changed without notice for quality improvement.