


ACCUMAX



PERFECTION REDEFINED

Accumax Heavy Metal Free Tubes

(15/50 ML Tubes)

Accumax Heavy Metal Free Tubes

(15/50 ml Tubes)

Heavy Metals Free centrifuge tubes are specifically manufactured to make sure that common trace metals that can interfere with experiments are below an ICP-MS detection level of one part per billion (<1PPB). Metals included in testing are: Boron, Aluminum, Phosphorous, Arsenic, Selenium, Cadmium, Lead, Chromium, Manganese, Iron, Nickel, Copper, Zinc, Antimony, Lithium, Silicon, Cobalt, Mercury, Magnesium, Potassium, Calcium, Sodium, Vanadium, Molybdenum



Product Features



Manufactured In Fully Automatic
Human Touch-free Class
100k Clean Room Facility With
Controlled Environment & Access

Withstands G-force
upto 20000



Produced From Swiss-made
Diamond Finish Mold For
Flawless 'Super Smooth'
Internal Surface

Certified below one part
per billion (<1PPB) trace
of heavy metals



DNase, RNase , Endotoxin ,
Pyrogen free

Tubes And Caps
Made With Virgin,
100% Medical-grade Resins



US FDA Approved
Virgin Polypropylene
Without Any Trace Of Metals

Tamper proof Graduation
& Large writing area



Ordering Information

Accumax Heavy Metal Free Tubes				
Cat. No.	Description	Packing Type	Quantity / Pack	Quantity / Case
ACT-15-B-T	15 ML Centrifuge Tube With Heavy Metal Free Flat Cap	Bulk	25 Tubes	25 Tubes X 20 packs
ACT-15-B-S-T	15 ML Centrifuge Tube With Heavy Metal Free Flat Cap, Sterile	Bulk	25 Tubes	25 Tubes X 20 packs
ACT-50-B-T	50 ML Centrifuge Tube With Heavy Metal Free Flat Cap	Bulk	25 Tubes	25 Tubes X 20 packs
ACT-50-B-S-T	50 ML Centrifuge Tube With Heavy Metal Free Flat Cap, Sterile	Bulk	25 Tubes	25 Tubes X 20 packs
ACT-50-K-B-T	50ml Self Standing Centrifuge Tube With Heavy Metal Free Flat Cap	Bulk	25 Tubes	25 Tubes X 20 packs
ACT-50-K-B-S-T	50ml Self Standing Centrifuge Tube With Heavy Metal Free Flat Cap , Sterile	Bulk	25 Tubes	25 Tubes X 20 packs
ACT-50-K-S-W-T	50ml Self Standing Centrifuge Tube With Heavy Metal Free Flat Cap, Sterile, Individually Wrapped	Individually Wrapped	25 Tubes	25 Tubes X 20 packs



**Caresource for Scientific and
Laboratory Equipment Trading CO LLC**

www.caresourceglobal.com
sales@caresourceglobal.com