

ACCUMAX



STURDY AND ATTESTED
MEDIA BOTTLES

Media Bottles
PETG/PET/PC

Media Bottle

Material: PETG/PET/PC

Volume: PETG : 30ml, 60ml, 125ml, 250ml, 500ml, 1L, 2L, 5L, 20L

PET : 125ml, 250ml, 500ml & 1L

PC : 30ml, 60ml, 125ml, 250ml, 500ml, 1L, 2L, 5L, 20L


Ergonomic Media bottles are suited for biopharmaceutical processing, storage of critical reagents, Media and various buffers used in laboratories. Manufactured in a clean room facility using high-quality laboratory-grade raw materials. Packaging is designed for clean room applications. The bottle and closure are secured with a tamper-resistant heat-shrink band, the bottle body is square shaped, making it convenient for storage purposes, The surface is smooth and shiny, easy to clean, and not easy to leave a residue. Our resin materials are chosen to minimize additives and potential leachables. We don't use any fillers or plasticizers.



Media Bottle

Secure, Dependable & Leakproof

When it comes to preserving valuable samples, compromise is not an option. Our range of media bottles is meticulously designed to offer unparalleled leakproof performance and deliver the utmost protection to samples



Sealed to Perfection

Threads on bottle and closures have continuous straight-shoulder semi-buttress threads ensuring a complete leakproof seal and preventing torque degradation



Product Features



Made of USP Class VI
Raw material, with low
leachable and extractable

Available in various sizes
ranging from 30ml to 20L



Manufactured in
Class 7 Clean room

Gamma sterile
product



Non Pyrogenic and
Non cytotoxic

Graduated for quick
user visibility



*Leakproof bottles &
HDPE/PP closures provided

PETG & PET can withstand
-70°C to 60°C ,
PC can withstand
121°C to -80°C



Heat-shrink band around closure
& neck for tamper-evident seal

* Leak proof testing is done using water, these tests using other liquids may not yield the same results.
Tested at 715mmhg pressure at room temp

Ordering Information

PETG			
Cat. No.	Product Description	Pack Size	Case Size
UAMB-30-PG-S	30ml Media Bottle with Cap, PETG, Sterile	24	96
UAMB-60-PG-S	60ml Media Bottle with Cap, PETG, sterile	24	96
UAMB-125-PG-S	125ml Media Bottle with Cap, PETG, sterile	24	96
UAMB-250-PG-S	250ml Media Bottle with Cap, PETG, sterile	24	48
UAMB-500-PG-S	500ml Media Bottle with Cap, PETG, sterile	24	48
UAMB-1000-PG-S	1L Media Bottle with Cap, PETG, sterile	12	24
UAMB-2000-PG-S	2 L Media Bottle with Cap, PETG, sterile	6	12
UAMB-30-PG-IW-S	30ml Media Bottle with Cap, PETG, Individually wrapped Sterile	1	46
UAMB-60-PG-IW-S	60ml Media Bottle with Cap, PETG, Individually wrapped sterile	1	40
UAMB-125-PG-IW-S	125ml Media Bottle with Cap, PETG,Individually wrapped sterile	1	40
UAMB-250-PG-IW-S	250ml Media Bottle with Cap, PETG, Individually wrapped sterile	1	54
UAMB-500-PG-IW-S	500ml Media Bottle with Cap, PETG, Individually wrapped sterile	1	70
UAMB-1000-PG-IW-S	1L Media Bottle with Cap, PETG, Individually wrapped sterile	1	50
UAMB-2000-PG-IW-S	2L Media Bottle with Cap, PETG Individually wrapped sterile	1	20
AMB-5L-PG-IW-S	5L Media Bottle with Cap, PETG , Individually wrapped, Sterile	1	4
AMB-5L-PG-IW	5L Media Bottle with Cap, PETG , Individually wrapped, Non Sterile	1	4
AMB-20L-PG-IW-S	20L Media Bottle with Cap, PETG , Individually wrapped, Sterile	1	1
AMB-20L-PG-IW	20L Media Bottle with Cap, PETG , Individually wrapped, Non Sterile	1	1
UAMB-30-PG-B-S	30ml Media Bottle with Cap, PETG, Bulk, Sterile	24	288
UAMB-60-PG-B-S	60ml Media Bottle with Cap, PETG, Bulk, Sterile	24	194
UAMB-125-PG-B-S	125ml Media Bottle with Cap, PETG, Bulk, Sterile	24	96
UAMB-250-PG-B-S	250ml Media Bottle with Cap, PETG, Bulk, Sterile	30	60
UAMB-500-PG-B-S	500ml Media Bottle with Cap, PETG, Bulk, Sterile	24	48
UAMB-30-PG-B	30ml Media Bottle with Cap, PETG, Non Sterile	24	288
UAMB-60-PG-B	60ml Media Bottle with Cap, PETG, Non Sterile	24	194

UAMB-125-PG-B	125ml Media Bottle with Cap, PETG, Non Sterile	24	96
UAMB-250-PG-B	250ml Media Bottle with Cap, PETG, Non sterile	30	60
UAMB-500-PG-B	500ml Media Bottle with Cap, PETG, Non sterile	24	48
UAMB-1000-PG-B	1L Media Bottle with Cap, PETG, Non Sterile	12	24
UAMB-2000-PG-B	2L Media Bottle with Cap, PETG, Non sterile	6	12
PET			
Cat. No.	Product Description	Pack Size	Case Size
AMB-125-PT	125ml Media Bottle with Cap, PET, Non sterile	24	96
AMB-250-PT	250ml Media Bottle with Cap, PET, Non sterile	30	60
AMB-500-PT	500ml Media Bottle with Cap, PET, Non sterile	24	48
AMB-1000-PT	1L Media Bottle with Cap, PET, Non sterile	12	24
AMB-125-PT-S	125ml Media Bottle with Cap, PET, Sterile	24	96
AMB-250-PT-S	250ml Media Bottle with Cap, PET, Sterile	30	60
AMB-500-PT-S	500ml Media Bottle with Cap, PET, Sterile	24	48
AMB-1000-PT-S	1L Media Bottle with Cap, PET, Sterile	12	24
AMB-125-PT-IW-S	125ml Media Bottle with Cap, PET, Individually Wrapped, Sterile	1	40
AMB-250-PT-IW-S	250ml Media Bottle with Cap, PET, Individually Wrapped, Sterile	1	54
AMB-500-PT-IW-S	500ml Media Bottle with Cap,PET,Individually Wrapped, Sterile	1	70
AMB-1000-PT-IW-S	1L Media Bottle with Cap, PET, Individually Wrapped, Sterile	1	50
PC MEDIA BOTTLE			
Cat. No.	Product Description	Pack Size	Case Size
UAMB-30-PC-S	30ml Media Bottle with Cap, PC, Sterile	24	96
UAMB-60-PC-S	60ml Media Bottle with Cap, PC, Sterile	24	96
UAMB-125-PC-S	125ml Media Bottle with Cap, PC, Sterile	24	96
UAMB-250-PC-S	250ml Media Bottle with Cap, PC, Sterile	24	96
UAMB-500-PC-S	500ml Media Bottle with Cap, PC, Sterile	24	48
UAMB-1000-PC-S	1L Media Bottle with Cap, PC, Sterile	12	24

UAMB-2000-PC-S	2L Media Bottle with Cap, PC, Sterile	6	12
UAMB-30-PC	30ml Media Bottle with Cap, PC, Non sterile	24	96
UAMB-60-PC	60ml Media Bottle with Cap, PC, Non sterile	24	96
UAMB-125-PC	125ml Media Bottle with Cap, PC, Non sterile	24	96
UAMB-250-PC	250ml Media Bottle with Cap, PC, Non sterile	24	96
UAMB-500-PC	500ml Media Bottle with Cap, PC, Non sterile	24	48
UAMB-1000-PC	1L Media Bottle with Cap, PC, Non sterile	12	24
UAMB-2000-PC	2L Media Bottle with Cap, PC, Non sterile	6	12
AMB-5L-PC-IW-S	5L Media Bottle with Cap, PC, Individually wrapped, Sterile	1	4
AMB-5L-PC-IW	5L Media Bottle with Cap, PC, Individually wrapped, Non Sterile	1	4
AMB-20L-PC-IW-S	20L Media Bottle with Cap, PC, Individually wrapped, Sterile	1	1
AMB-20L-PC-IW	20L Media Bottle with Cap, PC, Individually wrapped, Non Sterile	1	1

Packed in shrink-wrapped trays.



Caresource for Scientific and Laboratory Equipment Trading CO LLC

Tel: +971 4 385 6348

Mobile: +971 52 941 8683

Whatsapp: [Click Here](#)

E- mail: sales@caresourceglobal.com

Website: www.caresourceglobal.com

Al lotus bldg, 2nd Floor R-208, Salah Al Din St, Al khabaisi (main road), near abu baker metro station Deira – Dubai UAE