

A NEW DIMENSION TO ACCURACY


Accumax Serological Pipette

ACCUMAX



Manufactured From Virgin
Polystyrene In A Clean Room

Negative Graduation
Marks For Extra Capacity



Certified Free From
DNase, RNase, Cytotoxin
Free & Non-pyrogenic

Adhering To 10^{-6} SAL
(Sterility Assurance Level)
Quality Standards



Calibrated For Accurate
Dispensing Within
+/- 2% Volume Range

Filter Plug Prevents Both
Fluid And Liquid Vapors
Minimizing Contamination
Of Pipettor And Sample



Color Coding For
Easy Identification
Of Pipette Size



For more details visit us at www.caresourceglobal.com